

DIASONO•310

Modes

DIASONO·310

Modes

- ACNE
- Treatment Time : 18min

Frequency MHz	Strength	Time(min)	Treatment Period	# of Times	Treatment Target / Precautions
3/10	3MHz→30 10MHz→50	3	1~2 time(s)/week	8~10 times	<ul style="list-style-type: none">• Severe Acne• Redness caused by acne inflammation• Common Acne <p>※ Make sure to remove acne before treatment. ※ For acne treatments, you should explain in advance that client may experience more acne temporarily from 1 to 3 treatments.</p>
10	40	3			
3	60	3			
3/10	3MHz→30 10MHz→50	3			
10	40	3			
3	60	3			

DIASONO·310

Modes

- LIFTING
- Treatment Time : 18min

Frequency MHz	Strength	Time(min)	Treatment Period	# of Times	Treatment Target / Precautions
3/10	3MHz→30 10MHz→70	3	1~2 time(s)/week	8~10 times	<ul style="list-style-type: none">• Healthy skin that is not sensitive (Sensitive skin first do the REJUVENATION mode 3 times then use lifting mode)
10	80	2			
3/10	3MHz→50 10MHz→50	2			
3	60	2			
3/10	3MHz→30 10MHz→70	3			
10	80	2			
3/10	3MHz→50 10MHz→50	2			
3	60	2			

• Apply to Double chin/Lack of elasticity/Thick skin

• Increase elasticity by organizing facial lines

• Make sure to proceed with lymphatic drainage around the neck and ears before treatment.

※ If use high levels to sensitive skin, can cause skin troubles such as dead skin cells, itchiness, and redness.

Recommend the treatment after 3 adaptation periods with low level.

DIASONO·310

Modes

- ANTIAGING
- Treatment Time : 12min

Frequency MHz	Strength	Time(min)	Treatment Period	# of Times	Treatment Target / Precautions
3/10	3MHz→50 10MHz→50	3	1~2 time(s)/week	8~10 times	<ul style="list-style-type: none">• Young Skin Maintenance• Increased elasticity of sensitive skin• For sensitive skin, first do the REJUVENATION mode 1-2 times
10	50	3			
3/10	3MHz→50 10MHz→50	3			
10	50	3			

DIASONO·310

Modes

● REJUVENATION

● Treatment Time : 12min

Frequency MHz	Strength	Time(min)	Treatment Period	# of Times	Treatment Target / Precautions
3/10	3MHz→20 10MHz→40	2	1~2 time(s)/week	8~10 times	<ul style="list-style-type: none">• Activate skin regeneration• Improve skin texture• Improve skin tone• Pore Improvement• Increased moisturization of dry skin• Strengthen skin barrier easily reddening or prone to skin trouble• When skin regeneration, such as wound regeneration• Side effects from cosmetics, etc.• Skin with weak regenerative ability• Rapid regeneration of skin damage caused by treatment
10	50	2			
3/10	3MHz→30 10MHz→30	2			
3/10	3MHz→20 10MHz→40	2			
10	50	2			
3/10	3MHz→30 10MHz→30	2			

DIASONO·310

Modes

● WRINKLE

● Treatment Time : 12min

Frequency MHz	Strength	Time(min)	Treatment Period	# of Times	Treatment Target / Precautions
3/10	3MHz→30 10MHz→70	3	1~2 time(s)/week	8~10 times	<ul style="list-style-type: none">• Wrinkles around the eyes• Mouth-tailed wrinkles• Forehead wrinkles• Ear wrinkles• Thick wrinkles• Fine wrinkles <p>※ Slowly move 3mm per second as if ironing the treatment area.</p>
10	80	3			
3/10	3MHz→30 10MHz→70	3			
10	80	3			

DIASONO·310

Modes

- SEDATION
- Treatment Time : 12min

Frequency MHz	Strength	Time(min)	Treatment Period	# of Times	Treatment Target / Precautions
3/10	3MHz→40 10MHz→40	3	1~2 time(s)/week	8~10 times	<ul style="list-style-type: none">• natural redness• partially bounded redness
10	30	3			
3/10	3MHz→40 10MHz→40	3			
10	30	3			

DIASONO·310

Modes

- TONE UP
- Treatment Time : 12min

Frequency MHz	Strength	Time(min)	Treatment Period	# of Times	Treatment Target / Precautions
3/10	3MHz→30 10MHz→60	3	1~2 time(s)/week	8~10 times	<ul style="list-style-type: none">• Pigmentation due to acne Improvement• Skin tone Improvement• Expanded Pore Improvement• Scar Improvement
10	70	3			
3/10	3MHz→30 10MHz→60	3			
10	70	3			

DIASONO·310

Modes

- BRUISE&SWELLING
- Treatment Time : 9min

Frequency MHz	Strength	Time(min)	Treatment Period	# of Times	Treatment Target / Precautions
10	30	3			
3/10	3MHz→30 10MHz→30	3	-	1 time&more	<ul style="list-style-type: none">• Bruising/swelling areas after treatment• Bruising/swelling due to external impact• Treat application area for 5 seconds
10	50	3			<p>※ For swelling or bruising caused by filler treatment or fat transplantation. Do not apply.</p>

DIASONO•310

USER Mode

Guide

DIASONO·310

Modes

- USER Mode Guide

Treatment Target	Frequency	Strength	Time(min)	Treatment Period	# of Times
Laser Treatment(Fraxel, Thermage, etc) / thread lifting/ before needle treatment	3/10	3MHz→30 10MHz→40	12 min	-	1 time
Before Plastic Surgery	3/10	3MHz→30 10MHz→40	12 min	-	1 time
Ampoule (Solution,Booster) Penetration	3/10	3MHz→30 10MHz→40	12 min	-	1 time
Eczema	10	30	6 min	1~2time(s)/week	8~10 times
Atopic dermatitis	10	20	7 min	1~2time(s)/week	8~10 times
Keloid	10	50	7 min	1~2time(s)/week	8~10 times
Acne	3/10	3MHz→20 10MHz→40	12 min	1~2time(s)/week	8~10 times